


### Ordering Information

---



\* The Series 1000 Oxygen Monitor may be ordered with either RS-232C or RS-485.

- Configuration examples:
1. Series 1000 O<sub>2</sub> Monitor without options (*p/n 1000*)
  2. Series 1000 O<sub>2</sub> Monitor with remote sensor and RS-485 output (*p/n 1000-R-485*)
  3. Series 1000 O<sub>2</sub> Monitor with RS232C output and Battery Backup (*p/n 1000-232-B*)
  4. Replacement Sensor is 1SEN at a price of \$295

Note: The Calibration Fixture is recommended when the Series 1000 Oxygen Monitor is calibrated using a calibration gas.